## Advantage™

**Ultrasound Tissue Micking Phantom** 

Delivering Optimal Phantom Images

**Every Time** 

Long-lasting: Built with materials suitable for up to 10 years of continual use, the Advantage<sup>TM</sup> Phantom is meant for long-term application — with no loss in image quality

Reliable: Ultrasound users can easily verify the functionality of their equipment.

ACERTARA

**£**0

旕

Ø

53

ADVANTAGE

00000

HORIZONTAL

0

0

0

0

0

0

0

0

0

0

VERTICAL O O

0

0

AXIAL/LATERAL RESOLUTION

0

0

DEAD ZON

B C D 0 0 0 0

Gray Scale

HORIZONTAL

0 ± 10 mm

Lightweight: With a low physical weight of less than 5.5lbs/2.5kg, a small footprint, and a budget friendly price point, Advantage is the ideal QA phantom for the busy ultrasound lab.

**Ready to go:** Advantage requires no calibration or other downtime maintenance.

0

ACERTARA



## Acoustic Properties

- + Speed of sound: 1500±50 m/s
- + Attenuation: 0.6±0.15 dB/cm/MHz between 2 and 10 MHz (ambient temperature 21±3°C)

## Image Quality

+ Bright, easy-to-see targets with uniform homogenous backscatter field.

## Dual Scan

+ New dual scan surface design to easily accommodate linear, sector, endoscopic probes and transesophageal (TEE) probes.

See the Advantage<sup>™</sup> Phantom in action! Contact us today: 303.834.8413 | Sales@Acerteralabs.com